




A.




B.




C.




D.



E.



F.



The EVG CellProfiler pipeline workflow.

[A]. The original image (left) is masked using PathProfiler Tissue Segmentation Unet and used as input by CellProfiler 4.2.6; the graph (right) shows the tonal distribution in the digital whole-slide image on a RGB scale.

[B]. The input image is converted to a gray scaled image (left); the graph (right) shows the tonal distribution in the gray scaled image.

[C]. The gray scaled image is inverted, i.e. non-tissue will become black (left); the graph (right) shows the tonal distribution after inverting.

[D]. The tissue area is identified, as demarcated by the green line in the left image; the total tissue area size is calculated in pixels (right image) and tabulated (table).



[E]. The colors, i.e. stains, are unmixed using the original image (left): EVG (middle), and EVG_IDENTIFYING (HE, right). The graph (right) shows the tonal distribution of the EVG_IDENTIFYING image.

[F]. The EVG-positive objects, i.e. elastic tissue, are identified, white areas in the left image; the EVG-positive objects are demarcated by a green line in the middle image, areas that are excluded due to size (minimal size 1 pixels) are demarcated in magenta; the right image shows all the identified EVG-positive objects in random colors; the total number of identified objects is calculated and tabulated (table).

[G]. Finally, the data for each tile are saved in a comma-separated table, including meta-data such as tile positions, image location, object counts (there could be multiple patches of stained areas or tissue). The original image (top-left) is used to outline the EVG-positive objects. The tissue area (green), and EVG-positive objects (red) are all demarcated in the top-right image. The table (bottom-right) shows the areas occupied by each object class.

Sample used: AE9.T02-7170.EVG.TIF [Tile= X8000, Y24000]

G.



Objects or Image	Area Occupied	Perimeter	Total Area
EVG_Tissue	68181	33470.0	4000000
Tissue	1728224	16514.0	4000000